
Supplementary information

A metabolomics pipeline for the mechanistic interrogation of the gut microbiome

In the format provided by the
authors and unedited

Supplementary Methods

Mass spectrometry LC/MS methods

Instrumental and chromatographic settings.

Published C18 methods¹ and HILIC method² were implemented with minor modifications. The C18 positive method (ESI+) used mobile phase solvents (LC-MS grade) consisting of 0.1% formic acid (Fisher) in water (A) and 0.1% formic acid in methanol (B). The C18 negative method (ESI-) used mobile phase solvents consisting of 6.5 mM ammonium bicarbonate (Sigma) in water at pH 8 (A) and 6.5 mM ammonium bicarbonate in 95:5 v/v methanol:water (B). For both C18 methods, the gradient profile was from 0.5% B to 70% B in 4 minutes, from 70% B to 98% B in 0.5 minutes, and holding at 98% B for 0.9 minute before returning to 0.5% B in 0.2 minutes. The HILIC (Hydrophilic Interaction Liquid Chromatography) positive method (ESI+) used mobile phase solvents consisting of 0.125% formic acid and 10 mM ammonium formate (Sigma) in water at pH 3 (A) and 0.125% formic acid in 10 mM ammonium formate in 95:5 v/v acetonitrile:water (B). The gradient profile was held at 100% B for 2 minutes, from 100% B to 70% B in 5 minutes, holding at 70% B for 0.7 minute, from 70% B to 40% B for 1.3 minutes, holding at 40% B for 0.5 minutes, from 40% B to 30% B for 0.75 minutes, before returning to 100% B for 2.5 minutes and holding at 100% B for 4 minutes. The flow rate was 350 μ L per minute (C18 positive, C18 negative), and 400 μ L per minute (HILIC positive). The sample injection volume was 5 μ L (C18 positive, C18 negative), and 3 μ L (HILIC positive). LC separations were made at 40°C on separate columns fitted with a Vanguard pre-column of the same composition dedicated to each analytical method: Waters Acquity BEH 1.7 μ m particle size, 2.1 mm id x 100 mm length (C18), Waters Acquity BEH Amide 1.7 μ m particle size, 2.1 mm id x 150 mm length (HILIC). For all analytical methods, data were collected at a mass range of 70-1000 m/z at an acquisition rate of 2 spectra per second. Specific ion source parameters included Fragmentor (140V), Gas Temp (250°C), Sheath Gas Temp (200°C), and VCap (4000V).

Metabolomics sample preparation.

All samples were stored at -80°C until use, and were thawed on ice immediately before extraction. 200 μ L of serum and bacterial supernatant samples were used for extraction directly without dilution. Urine samples were diluted 1:20 in LC-MS grade water (Fisher) to reach a final volume of 200 μ L prior to extraction. Protein precipitation for bacterial, serum, and urine samples was conducted by adding 1 mL of extraction buffer (see composition below) in 100% methanol (LC-MS grade, Fisher) to 200 μ L of each sample in a 2 mL 96-well microplate (Fisher), sealed with a silicone mat (Agilent), and vortexed to mix. For feces and cecal contents, ~ 20 mg of feces or cecal contents were added to ~ 20 mg of acid-washed glass beads (150-212 μ m, Sigma) in a 2 mL autoclaved screw top vial. In the same vial, 600 μ L of water (LC-MS grade, Fisher) and 600 μ L of recovery buffer in 100% methanol (see composition below) were added. Fecal and cecal slurries were homogenized at 4°C using a Mini Beadbeater operating at 3,500 oscillations per minute for 5 minutes. For all sample types, samples were subsequently incubated at room temperature for 5 minutes, followed by centrifugation at 5,000 x g for 10 minutes. Two 440 μ L aliquots of the same supernatant were transferred to two separate 2 mL plates and dried under air in a Biotage TurboVap. One of these dried plates was sealed and archived at -80°C. The dried extracts were reconstituted in 200 μ L reconstitution buffer (see composition below) in 50% methanol in water (v/v, LC-MS grade, Fisher) by vortexing at max

speed for 5 seconds. Reconstituted sample extracts were centrifuged at 2,000 x g for 1 minute, and filtered through a 96-well Durapore PVDF 0.22-um filter plate (Millipore) into a 1 mL 96-well plate (Agilent) by centrifugation at 2,000 x g for 10 minutes. Plates were then sealed with 96-well cap mats (Agilent) and stored at -80°C until LC-MS analysis. QC samples were generated by pooling 4 µL from each well of the experiment into a single designated well on the same plate for LC-MS analysis.

The extraction buffer consisted of 4-Chloro-phenylalanine (6.8 µM, Carbosynth), Tridecanoic acid (6.8 µM, Sigma), and 2-Fluorophenylglycine (3.4 µM, SCBT) in 100% methanol. The reconstitution buffer included the internal standards: Phenylalanine-2,3,4,5,6-d5 (12.5 µM, CIL), Glucose-1,2,3,4,5,6,6-d7 (25 µM, CIL), Methionine-methyl-d3 (12.5 µM, CIL), 4-Hydroxyphenyl-d4-alanine (3.125 µM, CDN), Tryptophan-2,4,5,6,7-d5 (12.5 µM, CDN), Leucine-5,5,5-d3 (12.5 µM, CDN), N-Benzoyl-d5-glycine (6.25 µM, CDN), 4-Bromo-phenylalanine (12.5 µM, Sigma), Progesterone-d9 (3.125 µM, CIL), Di-N-octyl phthalate-3,4,5,6-d4 (12.5 µM, CDN), d19-Decanoic acid (12.5 µM, CDN), d15-Octanoic acid (25 µM, CDN), Indole-2,4,5,6,7-d5-3-acetic acid (12.5 µM, CDN), Carnitine-trimethyl-d9 (3.125 µM, CDN), and d27-Tetradecanoic acid (12.5 µM, CDN), in 50% methanol in water. The final concentration for each internal standard in these buffers was determined by choosing a concentration falling within its linear dynamic range as measured by each analytical method.

Data analysis

Custom bioinformatics: in vitro pipeline

The *in vitro* metabolomics pipeline was designed to enable comparison across experimental and analytical batches. Five main quality control steps were used (in order): removal of samples with high intra-replicate variability, removal of samples with substantial intra-experimental variation, internal standard and media blank selection and normalization, selection of measured molecular features (metabolites) on the basis of intra-replicate and intra-experiment variability, and removal of samples based on RF classification accuracy of taxonomy. Code is available in ‘in_vitro_pipeline.ipynb’, depicted in Extended Data Fig. 2, and summarized below.

Individual replicate groups were compared by Pearson’s correlation coefficient. Intra-replicate ion count (area under the curve) variability was generally low, with average correlation coefficient of 0.9-1.0 (log transformed data) depending on batch (Extended Data Fig 4b). Individual replicates were removed if their correlation with the other replicates was two standard deviations below the mean of all replicates. Multiple thresholds and removal strategies produced similar results. Entire batches (a set of samples that were harvested and ran on the instrument together) were then compared to identify outlier samples. Samples were compared via Principal Component Analysis (PCA). Samples distant from the centroid of the first two PCs were inspected at the chromatogram level, with focus on metabolites with high loadings on those PCs. In general, this strategy identified samples with substantial differences in single metabolites due to technical artifacts (e.g. poor peak integration, non-detection, column pressure fluctuation, etc.).

LC-MS measurements are frequently confounded by matrix effects (differences in ionization efficiency for the same metabolite present in different biological sample types) and changes in instrumental sensitivity over time. To minimize these effects, we normalized samples based on a shared set of internal standards (IS). After sample removal, IS ion counts were compared within and between experiments on a per-analytical method basis. Based on matrix-

effect experiments, we assumed that highly variable ($>10X$ difference in count between samples, intra-experiment) IS were not representative of ionization efficiency and were excluded from analysis. The remaining panel of IS (2-Fluorophenylglycine, 4-Bromo-phenylalanine, 4-Chloro-phenylalanine, d19-Decanoic acid, Methionine-methyl-d3, N-Benzoyl-d5-glycine, Indole-2,4,5,6,7-d5-3-acetic acid, Phenylalanine-2,3,4,5,6-d5, and Tridecanoic acid) were summed (per sample) and used to compute a multiplicative scaling factor to equalize each sample ($\frac{\sum IS}{\sum_{sample} IS}$). This was done on a per-analytical method basis. For example, the correction factor for C18 positive method was computed from only IS detected in C18 positive method and applied to only detected features using this method. Counts of all metabolites (per sample and per method) were multiplied by the samples (method-specific) scaling factor prior to further analysis. We explored other normalization strategies including a scaling factor based on a weighted sum of IS, but chose total IS sum because it minimized intra-replicate coefficients of variation (CV).

In addition to IS normalization, we computed several other data transformations with the goal of allowing inter-experiment comparisons (Supplemental Table 7). IS-based correction helped to normalize average molecular feature intensity across runs but could not correct for all inter-experiment variation or biological differences, e.g., differences in molecular feature count due to media batch. For *in vitro* data used in the figures (e.g., Fig. 3) and in the Metabolomics Data Explorer, the data were subjected to the following transformation to mitigate these concerns. First, a ‘delta’ matrix was constructed (IS-corrected supernatant samples minus IS-corrected media blanks on a per-molecular feature basis). For each replicate group (i.e. three biological replicates), molecular feature counts were eliminated if the replicate group average was less than 3000 counts or the delta was less than 3000 counts (absolute value). The purpose of this filter was to eliminate peak counts that were likely artifactual (e.g., due to integration of a non-peak or a peak with low signal-to-noise ratio). This filter eliminated 9.9% of measured values. Different filters (1000-10000 count) were explored with minimal impact on the data since almost all of these counts were random noise or very low quality peaks. Second, molecular features that varied by more than 10X (by ion count) within a replicate group were eliminated. This threshold was chosen as the elbow of the distribution of (per feature) per replicate group maximum count divided by minimum count. This eliminated 1.0% of the measured values. Finally, we used this filtered data and computed a fold change matrix. Values in the fold change matrix are the count of a molecular feature in a supernatant sample divided by the count of that feature in the blank media that was used to grow that supernatant. This final matrix enables inter-experimental comparisons while minimizing technical and biological batch effects. All intermediate transformations and the raw untransformed data are available in Supplementary Table 7.

Because many metabolites were detected by more than one analytical method, we next identified a strategy for selecting which method(s) to prefer for fold-change data on a per-metabolite basis. In this step, columns of the data matrix were transformed from detected molecular features detected in multiple methods (redundant) to a single metabolites (non-redundant). For each molecular feature detecting the same metabolite we computed the CV (calculated on a per-replicate basis) for all replicates with those features. We compared the distributions of CVs and the linear correlation between molecular features. In cases where correlation was low, (Pearson $r < \sim 0.4$) we selected the feature with the better CV distribution (smaller mean, positively skewed) and ignored the other molecular feature. In cases where correlation was high, we averaged the fold changes. When a metabolite was detected by only one

molecular feature it was included regardless of how well that feature was correlated with the other detecting features.

The final quality control step evaluated the similarity of supernatant samples to all other samples from the same phylum. Preliminary work with our data showed that phylum level differences were evident in our supernatants, and they could be easily identified by random forest classifiers (RFs). Using a test:training split of 1:2, RF classifiers were trained to predict the phylum of the microbe generating each supernatant sample (see ‘sample_qc_classification.ipynb’). The phylum of each test sample was then predicted and misclassifications inspected. This analysis was repeated with different parameters (e.g., binary classification or multiclass classification) and with different input data (e.g., fold-change or ion count data). Consistently, misclassified samples were individually inspected for phyla-specific features (e.g., samples were contaminated by Indolepropionic acid-producing *C. sporogenes* from Phylum Firmicutes). Samples with significant counts (>10,000 count average across all replicates) of phyla-specific metabolites (from a different phyla) that did not replicate on repeated growth were eliminated.

Here we have presented this methodology as linear, but it should be noted that the filtration process was iterative with later steps informing earlier filtration steps. As an example, the PCA and RF analysis steps revealed molecular features that were inconsistent across runs and had to be eliminated in the initial data curation. In total, this pipeline represents a conservative approach more likely to eliminate ‘good’ samples than allow samples containing substantial artifacts (biological or technical) to remain in the data. Cumulatively this pipeline eliminated 8.13% of samples and 10.9% of measured counts in the retained samples.

Custom bioinformatics: in vivo pipeline.

The *in vivo* pipeline integrated data across mono-colonization, community, and conventional mouse experiments and three analytical methods to generate a unified, metabolite fold-change matrix. First, the sample database and MS-DIAL output files from the experiments were parsed and integrated into a single data matrix of raw ion counts across the three analytic methods. Based on the reference library, a small subset of metabolites exhibited dual peaks specific to an analytical method. For these metabolites, the raw ion counts from these two peaks were summed in the raw ion count matrix. Next, a shared set of internal standards were used to normalize across experiments to account for inter-experimental variations in instrument sensitivity. As the qTOF instrument was in a shared core facility that ran near full capacity, not a dedicated instrument for these metabolomics methods, run-to-run variation in sensitivity was an important issue. To account for these variations, the raw ion count from each metabolite was normalized to the sum of internal standard ion counts specific to each experiment (e.g., community experiment) and sample type (e.g., serum). The fold change matrix was generated next by calculating the relative fold change between metabolite ion counts detected in colonized mice vs. germ-free controls for each experiment and sample type. In particular, cecal samples from the community experiment were processed differently from other experiments, because the cecum from each mouse was trisected and three cecal samples were harvested per mouse. In this case, before calculating fold-change values, the three cecal samples for each mouse were collapsed into a single sample row in the ion count matrix by averaging the ion counts for each metabolite. Lastly, when a metabolite was detected by multiple analytical methods, its fold-change values were averaged among the preferred detection methods, which were determined based on the consistency (CVs) of detection among replicates for each method (see *in vitro* pipeline).

The final fold-change matrix was used for several analyses: 1) PCA analyses to separate metabolomic profiles by sample types or colonization states, 2) statistical calculation to identify significantly regulated metabolites (≥ 4 fold) using *t*-test with Benjamin-Hochberg corrections for multiple comparisons ($P < 0.05$), and 3) violin plots to visualize metabolites associated with a specific colonization state in a given sample type. All steps detailed here are explained in depth in the Jupyter notebook “mouse_data_analysis.ipynb”. The final fold-change matrix was used to construct Fig. 4 and Extended Data Figs. 8e, 9, 10.

Fold-change matrices generated by both *in vitro* and *in vivo* pipelines can be interactively accessed in the web-based Metabolomics Data Explorer (https://sonnenburglab.github.io/Metabolomics_Data_Explorer)

Bacterial sequencing and phylogenetics

Purity analysis

Culture purity was assessed with 16S amplicon sequencing of the V4 (515f, 806r) region using Golay barcoded primers following the protocols of the Earth Microbiome Project (<https://earthmicrobiome.org/>). Briefly, DNA from bacterial cultures was extracted using DNeasy PowerSoil HTP 96 Kit (Qiagen). PCR conditions were as follows: initial denaturation 3 minutes at 94°C, followed by 35 cycles of denaturation (45 seconds at 94°C), annealing (60 seconds at 50°C), extension (90 seconds at 72°C), followed by a final elongation step of 10 minutes at 72°C. PCR products were cleaned using UltraClean 96 PCR Cleanup Kit (Qiagen) and quantified with the Quant-iT PicoGreen dsDNA Assay Kit (Invitrogen). Cleaned products were pooled in equimolar amounts and sequenced using 300-bp paired end reads on a MiSeq (Illumina). Illumina reads were demultiplexed and sequence variants were determined using QIIME 1.9.1³ and DADA2⁴. Sequence variants in each sample were compared against the expected sequence of the grown strain taken from NCBI. A sequence match was considered to be at least 251 nucleotides (nt) out of 253 nts from the NCBI reference. If more than 5/1000 reads did not match, the culture was considered contaminated. Some exceptions for nt match thresholds were made for strains with low quality 16S sequences in the NCBI database (e.g. *Tyzzarella nexilis*).

Phylogenetics

Phylogenies were constructed with both V4 amplicons generated by Illumina sequencing (above) and with near full-length 16S sequences identified via automated and manual search of NCBI. 16S sequences >1200 nts in length were considered for full-length phylogeny reconstruction and strains were excluded if a full-length sequence could not be found. In both cases, phylogeny was constructed using Clustal Omega with default parameters using the EMBL-EBI web interface⁵.

Because strains frequently had more than one associated V4 16S variant, we created a modified phylogenetic tree with a one-to-one correspondence between tips and strains used in this study that was used in the distance calculations described below (see ‘phylogenetic_tree_modifications.ipynb’). Trees were visualized with iTOL⁶.

Distance comparisons and classifiers

Distance calculations

To generate distance matrices used for Fig. 2a, and Extended Data Fig. 7a, d-i, we followed the following procedure (detailed in ‘distance_comparisons.ipynb’). First, we averaged the measured values for metabolites for each strain across all experiments and log transformed the data. For each pair of strains, we calculated the Euclidean distance between the metabolite vector of those strains, excluding metabolites that were not measured in both (e.g., metabolites with NaNs in one or both vectors were removed prior to calculation). We also computed the tip-to-tip distances for all strains using the V4 region and full-length phylogenetic trees using scikit-bio. We also computed metabolomic distance matrices using input metabolomic data with various transformations and levels of filtration. Specifically, we used count data, unfiltered fold change data, and fold change data with media blank, delta, and variance ratio filtration. The outcome of the following analyses was unaltered with these different input data sets suggesting our results are robust and not the product of a well-chosen transformation.

The heatmap in Extended Data Fig. 6a with 158 mega-medium grown strains hierarchically clustered based on their metabolomic distance was generated as described in ‘in_vitro_heatmap_scatter_plot.ipynb’. The *in vitro* metabolite fold-change matrix (log2 transformed) was used. On the y-axis, taxonomies, or bacterial strains, are clustered by their metabolomic distance. We first computed the metabolomic distance matrix using the Euclidean distance between fold-change values as the distance metric. We then hierarchically-clustered (Ward’s method) the taxonomies and constructed the heatmap according to the resulting taxonomy order. On the x-axis, metabolites are also hierarchically-clustered (Ward’s method) using the Euclidean distance between the fold-change (log2 transformed) values across all taxonomies.

Generation of metabolomic trees weighted based on chemical similarity between the metabolites in the *in vitro* dataset (Extended Data Fig. 7b, c) used the procedure as described in ‘invitro_metabolite_similarity.ipynb’. We first computed a chemical distance matrix using pairwise Tanimoto 2D structural similarity scores (PubChem database scores range from 0 to 100 as the most similar) among all uniquely detected (non-coeluting) metabolites in the *in vitro* dataset. Based on this chemical distance matrix, we performed hierarchical clustering on the metabolites and assigned weights to each metabolite such that the sum of weights within each cluster of closely related metabolites equal 1. Lastly, we constructed weighted metabolomic trees using metabolomic profiles of bacterial taxa that are hierarchically-clustered based on their weighted Euclidean metabolomic distance.

To compare the weighted and unweighted metabolomic trees, we performed the Mantel test to correlate the two distance matrices computed on the weighted and unweighted Euclidean distance. We found that the matrices were highly correlated despite the addition of weighting. To transfer the metabolomic trees to iTOL for visualization, we converted the metabolomic clustering (expressed as a linkage matrix) into Newick trees, which iTOL supports as input. Supporting code is described in ‘invitro_metabolite_similarity.ipynb’.

Correlation between labels in phylogenetic and metabolomic tree

To understand if the relative position of the labels (strains and associated taxonomies) were similar between the phylogenetic and metabolomic tree we used the following procedure (Extended Data Fig 7a). For each tree we calculated the 4th percentile of the tip-to-tip distance (TTD) matrix. We designated this value as the ‘radius’ of a ‘neighborhood’ such that for a given tip, any other tip with pairwise distance less than the radius was considered part of the given tips neighborhood. For each tip, we compared the neighborhoods of that tip in the phylogenetic and

metabolomic tree. We calculated the overlap between these neighborhood sets at each taxonomic level. We used the number of tips in the phylogenetic neighborhood as the denominator to calculate the fraction of the labels seen in the metabolomic neighborhood that were found in the phylogenetic neighborhood. For example, if the phylogenetic neighborhood contained 5 Bacteroidetes and the metabolomic neighborhood 10 Bacteroidetes, the fractional representation would be 1.0 at the phylum level. If the phylogenetic neighborhood contained *B. ovatus* strain 1, *B. ovatus* strain 2, and *B. xylanisolvens* strain 1 and the metabolomic neighborhood contained *B. ovatus* strain 1, *B. ovatus* strain 10, and *P. distasonis* strain 1, the overlap would be 1/3 at the strain level and 2/3 at the species level. We repeated this procedure for all tips in the tree, and for definitions of neighborhood radius ranging from the 1st to 10th percentile of distances (in 1% increments). The results shown in Extended Data Fig. 7a are representative and are from a radius equal to the 4th percentile in each TTD respectively.

To calculate how significant the observed fractional coverage was, we used a permutative approach. For each tip, we calculated the true fractional overlap between neighborhoods. We then shuffled the labels of the TTD 1000 times, and recorded the number of times a fractional overlap larger than the observed occurred. The p-value (shuffled fractional overlap > observed fractional overlap / 1000) is reported in Extended Data Fig. 7a. The code for this section can be found in the github repository under “distance_comparisons.ipynb”.

Distance comparisons

For Extended Data Fig. 7a, d-i we compared the phylogenetic and metabolomic distances. Patterns were consistent across multiple different transformations of the data including data type (fold change or count), data filtration (no filter or variance and intensity filter), and NaN removal (no removal or replacement by 1). The LOESS regressions of the data matrices were done using *scipy*.

Random forest classification

To establish if there were consistent metabolomic signals of taxonomic labels (e.g., Bacteroidetes or Firmicutes) we used random forest (RF) classifiers. For Fig. 3a, 25 RFs were trained on each taxonomic level (e.g. to predict phylum labels) using metabolomic data (see ‘figure_3_classifiers.ipynb’). Fold change data were used to avoid batch effects of raw ion counts, though classification accuracy was not significantly different using count data. Data were split 1:2 test:training without class balancing due to the over-representation of Bacteroidetes and Firmicutes in the data.

For identification of Bacteroidetes specific metabolites, we trained RFs to predict phylum level labels of Bacteroidetes and Other (binary classification, Extended Data Fig. 8a, top panel). This strategy identified Bacteroidetes specific metabolic behaviors with high accuracy leading to ~97% classification accuracy on average. We built 50 forests and averaged the feature importance from all forests to produce the overall importance data in Extended Data Fig. 8a, bottom panel.

To avoid overfitting both of the above tasks, we significantly constrained the number of total decisions a tree could make (maximum depth ≤ 5) and eliminated features that were present in less than 30% of samples. Other thresholds yielded similar results, though below this depth, trees lost accuracy.

Mega Medium (MM) preparation protocol

This recipe was adapted from a previous publication⁷ and from a personal communication of N. McNulty. The starch was used as indicated in the supplemental methods.

Steps

1. Add dry ingredients, liquid ingredients, and mix well (table 1).
2. pH to ~7.0. In practice, this takes slightly less than 2.5 mL 10M KOH per 500 mL media. After pH measurement, add $3.6 - X$ mL water where X is mL KOH added.
3. If plates are needed, add 7.5 g agar per 500 mL of media.
4. Autoclave cycle 4 (250 F for 25 minutes at 20 PSI, 20 minutes dry).
5. As soon as autoclaving starts, remove vitamin mix from freezer to let thaw. Once thawed, make vitamin solution and filter sterilize (table 2).
6. After media has cooled, sterilely add vitamin solution.

Table 1

Component	Amount (in 500 mL)	Final conc.	Mfr.	Vendor (cat. #)
Milli-Q water (dH ₂ O)	410 mL			
1 M Potassium phosphate buffer, pH 7.2	50 mL	10% (v/v)		
TYG salts solution	20 mL	4% (v/v)		
Tween 80	1 mL of 25% (v/v)	0.05% (v/v)		
SCFA supplement	1.4 mL	0.28% (v/v)		
0.8% (w/v) CaCl ₂	500 µL			
FeSO ₄ ·7H ₂ O	500 µL of 0.4 mg/mL			
Resazurin	2 mL of 0.25 mg/mL	0.000001% (w/v)	Sigma	Sigma (R2127)
Trypticase Peptone (BBL)	5 g	1% (w/v)	BD	BD (211921)
Yeast Extract (Bacto)	2.5 g	0.5% (w/v)	BD	BD (212750)
Meat extract	2.5 g	0.5% (w/v)	Sigma	Sigma (70164)
L-Cysteine hydrochloride	0.25 g	0.05% (w/v)	Sigma	Sigma (C1276)

D-(+)-Glucose	1 g	0.2% (w/v)	Sigma	Sigma (G8270)
D-(+)-Cellobiose	0.5 g	0.1% (w/v)	Sigma	Sigma (C7252)
D-(+)-Maltose monohydrate	0.5 g	0.1% (w/v)	Sigma	Sigma (M5885)
D-(-)-Fructose	0.5 g	0.1% (w/v)	Sigma	Sigma (F0127)
Soluble starch	0.25g	0.05% (w/v)		

Table 2

Component	Amount(in 500 mL)	Final conc.	Mfr.	Vendor (cat. #)
Vitamin K solution	500 μ L of 1 mg/mL	0.000001% (w/v)	Sigma	Sigma (M5625)
Trace Mineral Supplement	5 mL	1% (v/v)	ATCC	ATCC (MD-TMS)
Vitamin Supplement	5 mL	1% (v/v)	ATCC	ATCC (MD-VS)
Histidine-Hematin	500 μ L			

STOCK SOLUTION RECIPES

1 M potassium phosphate buffer, pH 7.2

1. Prepare 1 M KH_2PO_4 (monobasic).
68.045 g KH_2PO_4 (anhydrous, f.w. = 136.09) in Milli-Q water to 500 mL
2. Prepare 1 M K_2HPO_4 (dibasic).
174.18 g K_2HPO_4 (anhydrous, f.w. = 174.18) in Milli-Q water to 1 L
3. Add monobasic to dibasic to achieve pH 7.2.
(You typically need ~430 mL monobasic added to 1 L dibasic)

Tip: When preparing this solution, add the potassium phosphate powders to an actively stirring volume of water. Attempting to add water to a large mass of powder will result in the formation of a difficult to dissolve clump of material at the bottom of the bottle.

Vitamin K solution

Dissolve 40 mg menadione (Vitamin K₃, Sigma M5625) in 40 mL 100% EtOH.

TYG salts solution

$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (Sigma 230391) 0.5 g
 NaHCO_3 (Sigma S5761) 10.0 g
 NaCl (Sigma S7653) 2.0 g

Milli-Q water to 1 L

FeSO₄·7H₂O (0.4 mg/mL)

Dissolve 40 mg FeSO₄·7H₂O (Sigma F8633) in 100 mL Milli-Q water.

0.8% (w/v) CaCl₂

Dissolve 0.4 g CaCl₂·2H₂O (Sigma C7902) in 50 mL Milli-Q water.

Resazurin anaerobic indicator (0.25 mg/mL)

1. Dissolve 25 mg resazurin (Sigma R2127) in 100 mL distilled H₂O.
2. Store protected from light at 4 °C.

Histidine-Hematin

1. Prepare 0.2 M histidine, pH 8.0
 - Mix 4.2 g Histidine-HCl monohydrate (Sigma H7875) in 80 mL Milli-Q water.
 - Adjust the pH from 4 to 8 with 10 N NaOH (the histidine will go into solution as the pH rises).
 - Bring the final volume to 100 mL with Milli-Q water.
2. Mix 12 mg hematin (Sigma H3281) with 10 mL of 0.2 M histidine, pH 8.0. Dissolve by end-over-end rotation or vigorous shaking for several hours. Filter-sterilize using 0.2 µm filter.

SCFA supplement

Acetic acid, glacial (Sigma A6283)	17 mL
Propionic acid (Sigma P5561)	6 mL
Butyric acid (Sigma B103500)	4 mL
Isovaleric acid (Sigma 129542)	1 mL

Salyer's Minimal Medium (SMM) preparation protocol

This recipe was adapted from multiple references^{8,9} that give slightly different concentrations for various components. The media was prepared and used as described in the supplemental methods.

Mineral salts (Table 3 below) and hematin solution can be prepared, filter sterilized, and stored at 4°C for at least 12 months. Carbon sources and nitrogen sources (critically if using glutamine) should be prepared the day of use.

Steps:

1. Prepare mineral salts, hematin, vitamin K3 solution, iron sulfate solution, and vitamin B₁₂ solution.
2. Prepare SMM base according to Table 1 below. Filter sterilize and place into anaerobic chamber for reduction. May take 48 hours given paucity of reductants.
3. Prepare carbon, nitrogen, and sulfur sources immediately prior to use, filter sterilize, and add to reduced SMM base according to Table 2 below.

Table 1 – SMM Base

Component	Amount (in 100 mL)
dH ₂ O	60.84 ml
1 M KPO ₄ , pH 7.2	10 ml
Vitamin K3 solution, 1mg/ml	0.01 ml
FeSO ₄ ·7H ₂ O, 0.4 mg/ml	1 ml
Vitamin B ₁₂ , 0.01 mg/ml	0.05 ml
Na ₂ CO ₃	0.1 g

Table 2 – SMM Complete

Component	Amount (in 100 mL)
SMM base	71.9 ml
Carbon source (50x)	2 ml
Mineral salts	5 ml
Reduced sulfur source (100X)	1 ml
Hematin solution	0.1 ml
Nitrogen source (5X)	20 ml

Table 3 – Mineral salts

Component	Amount (in 1L)
NaCl	18 g
CaCl ₂ .2H ₂ O	0.53 g
MgCl ₂ .6H ₂ O	0.40 g
MnCl ₂ .4H ₂ O	0.20 g
CoCl ₂ .6H ₂ O	0.20 g
dH ₂ O	1 L

Carbon source at a final concentration of 27.7 mM (50X=1.385M). Here we used Glucose

Component	Amount (in 50mL)
dH ₂ O	50 ml
Glucose	12.48 g

Hematin at a 100X concentration. To prepare this, first mix 0.4g NaOH in 100ml dH₂O, then add 0.5g Hemin or Hematin.

Component	Amount (in 100ml)
0.1 M NaOH (0.4g NaOH in 100ml dH ₂ O)	100 ml
Hemin	0.5 g

Reduced sulfur source at a final concentration of 4.12 mM (100X = 412.7mM). Here we used both cysteine and sodium sulfide.

Component	Amount (in 20mL)
dH ₂ O	20 ml
Cysteine	1 g
dH ₂ O	20 ml
Na ₂ S	0.644 g

Nitrogen sources at a final concentration of 10mM (5X = 50 mM). Ammonium sulfate was assumed to produce two equivalents and all others one equivalent of NH_4 .

Component	Amount (in 30mL)
dH ₂ O	30 ml
(NH ₄) ₂ SO ₄	0.1 g
dH ₂ O	30 ml
Glutamate	0.22 g
dH ₂ O	30 ml
Glutamine	0.22 g
dH ₂ O	30 ml
Asparagine	0.2 g

References

1. Wikoff, W. R. *et al.* Metabolomics analysis reveals large effects of gut microflora on mammalian blood metabolites. *Proc. Natl. Acad. Sci.* **106**, 3698–3703 (2009).
2. Showalter, M. R. *et al.* Obesogenic diets alter metabolism in mice. *PLoS One* **13**, 1–17 (2018).
3. Caporaso, J. G. *et al.* QIIME allows analysis of high-throughput community sequencing data. *Nat. Methods* **7**, 335–336 (2010).
4. Callahan, B. J. *et al.* DADA2: High-resolution sample inference from Illumina amplicon data. *Nat. Methods* **13**, 581–583 (2016).
5. Sievers, F. *et al.* Fast, scalable generation of high-quality protein multiple sequence alignments using Clustal Omega. *Mol. Syst. Biol.* **7**, (2011).
6. Letunic, I. & Bork, P. Interactive Tree of Life (iTOL) v4: Recent updates and new developments. *Nucleic Acids Res.* **47**, 256–259 (2019).
7. Wu, M. *et al.* Genetic determinants of in vivo fitness and diet responsiveness in multiple human gut Bacteroides. *Science* **350**, (2015).
8. Shepherd, E. S., Deloache, W. C., Pruss, K. M., Whitaker, W. R. & Sonnenburg, J. L. An exclusive metabolic niche enables strain engraftment in the gut microbiota. *Nature* **557**, 434–438 (2018).
9. Varel, V. H. & Bryant, M. P. Nutritional features of Bacteroides fragilis subsp. fragilis. *Appl. Microbiol.* **28**, 251–7 (1974).

Supplementary Table Guide

Supplementary Table 1

Title: Mass spectrometry compound m/z-RT reference library

Table columns:

comp. info columns	dname, compound, peak, pubchem_CID, mz (precursor accurate mass/charge), rt (retention time identified on the qTOF instrument), adduct, molecular_formula, monoisotopic mass, QE_rt (retention time identified on the Q Exactive instrument), QE_ms2 (“*” denotes ms2 spectra generated on QE, spectra reported in Supplementary Table 3), qTOF_ms2 (“*” denotes ms2 spectra generated on qTOF, spectra reported in Supplementary Table 3)
comp. class. columns	kingdom, superclass, class, subclass
comp. id columns	canonical_smiles, inchikey
analytical_method	“mode”: C18 positive, C18 negative, HILIC positive

Supplementary Table 2

Additional info:

MS/MS validation of biological metabolites and library standards, matrix effects on ion count and RT, and linear range of individual metabolites

Worksheet: ms2_validation

comp. info columns	library_dname, library_compound, library_mz, library_rt, library_adducts
best_matched_ce	biological compound or library standard, the collision energy (10, 20, and/or 40 eV) that produces the best-matched spectra to the MoNA database based on similarity scores
matched_compound	biological compound or library standard, the compound name of the matched spectrum from the MoNA database
ms/ms_similarity	biological compound or library standard, the similarity score computed by comparing MS/MS spectral fragments’ m/z and intensity between our spectra and the MoNA spectra, using a build-in spectral comparison algorithm in MoNA
analytical_method	C18 positive, C18 negative, HILIC positive
sample_type	applies to biological compound columns only, includes bacterial_supernatant, fecal_cecal (fecal and cecal pooled), serum, and urine

Worksheets: matrix_effects_retention_time, matrix_effects_ion_count, linear_range

comp. info columns	library_dname (or dname), library_compound (or compound), library_mz, library_rt, library_adduct
analytical_method	or “mode”: C18 positive, C18 negative, HILIC positive
sample cond. columns	Matrix effects-specific columns: Control, Control_spikein, Mouse_feces, Mouse_feces_spikein, Mouse_urine, Mouse_urine_spikein, Mouse_serum, Mouse_serum_spikein, Human_serum, Human_serum_spikein, Bacterial_media, Bacterial_media_spikein
conc. columns	linear range-specific columns: numerical values from 0.001 to 200 (uM)
lower_bound	the lower concentration at which the linear range is calculated, concentration indices (0 = 0.0011 and 11=200) are used for range annotation
upper_bound	the upper concentration at which linear range is calculated, concentration indices (0 = 0.0011 uM and 11 = 200 uM) are used for range annotation

m	slope of the linear regression line ($y = 10^b x^m$)
b	y-intercept of the linear regression line
r ²	coefficient of determination for the linear regression line

Supplementary Table 3

Title: MS/MS spectra library constructed on qTOF and QE instruments

Additional info:

Spectra library file containing MS/MS reference spectra from authentic standards collected on two distinct instruments (750 compounds, Agilent qTOF 6545 and 773 compounds, Thermo Orbitrap Q Exactive HF). Three standard collision energies (10, 20, and 40 eV) were collected for each compound on the qTOF instrument, and one normalized collision energy (20-30-40% NCE) was collected on the QE instrument.

Comp. info rows	Name, Precursor_mz, Precursor_type, Spectrum_type, InChIKey, SMILES, Formula
Instrument info rows	Instrument_type, Instrument, Ion_mode, Collision_energy
Spectra info rows	Num Peaks, m/z and intensity pairs separated by ‘,’

Supplementary Table 4

Title: Inter-instrumental retention time shift and correction

Additional info:

Calculation of retention time (RT) difference and non-linear RT correction between two mass spectrometry instruments: library instrument qTOF 6454 vs. a second instrument qTOF 6530, or library instrument qTOF 6454 vs. a second type of instrument orbitrap Q Exactive (QE), for a shared panel of reference library standards spanning diverse RTs.

Worksheets: qTOF_RT_correction and QE_RT_correction

Table columns:

comp info columns	library_dname, library_compound, library_mz, library_rt, library_adducts
corrected_rt	library rt corrected by polynomial transformation of the library based on inter-instrumental RT shifts of 10-20 robustly detected metabolites (e.g. internal standards)
measured_rt	RT measured on the second instrument
measured vs. library	measured RT on the second instrument minus library RT
measured vs. corrected	measured RT on the second instrument minus corrected library RT
analytical_method	C18 positive, C18 negative, HILIC positive

Supplementary Table 5

Title: *in vitro* sample metadata

Additional info: this table contains the metadata associating each bacterial supernatant sample to the culture and LC-MS conditions that produced the metabolomic data for that sample. It also contains data on the optical density (OD) measurements.

Worksheet: mf

sample_id	stable sample ID associating a single supernatant or media blank measurement with its metabolite measurements.
-----------	--

sample_type	supernatant, media_blank, qc, istd_blank, blank; the different types of sample run on the LC-MS instrument
clean_by_16s	records that the supernatant is from a pure and verified stock
culture_source	the culture tube used
matrix_tube	the matrix tube used
taxonomy	the genus and species of microbe cultured (if a supernatant sample)
experiment	identifier for growth batch
od_sample_id	the unique sample ID associating the sample to its replicate groups OD measurement in the 'od' worksheet
storage_loc_plate	what plate a sample is stored in
storage_loc_well	what well a sample was stored in
preculture_date	date of supernatant preculture
preculture_time	length of preculture in hours (approximate)
subculture_time	length of subculture in hours
media	media used; for abbreviations see supplemental table 6 worksheet 'media'
mode	LC/MS analytical method: c18positive, c18negative, or hilicpositive
chromatography	c18 or hilic
ionization	positive or negative
lcms_run_date	date samples were run
run_designation	batch within an LC-MS run
run_order	order of samples in run

Worksheet: ods

sample_id	unique sample ID associating the optical density measurements to a specific sample
od_time_step	time (hours) between each OD reading
m0001...	all these headers correspond to sequential OD measurements
taxonomy, experiment, preculture_date, preculture_time, subculture_time, media	are the same as the 'mf' worksheet

Worksheet: sample_link

c18positive, c18negative, hilicpositive repeats	sample IDs that measure the same physical sample repeated technical injections
---	--

Supplementary Table 6

Title: culture and strain information

Additional info: this table contains all the information on the taxonomy, NCBI accession numbers, 16S sequences and the like on the bacterial strains used in this manuscript.

Worksheet: strains_with_data

This worksheet lists the bacterial strains that correspond to data in Supplementary Table 7

Worksheet: cultures

sonnenburg_culture_num	a unique identifier for a glycerol stock tube, shared with Supplementary Table 5
taxonomy	taxonomy of the glycerol stock
from_matrix_tube	matrix tube used for particular cultures, shared with Supplementary Table 5
tube_top_nums	matrix tube identifier
clean_by_16s	if culture is clean by 16S analysis

16s_feature	the 16S v4 region feature name (or names) for this strain, used for constructing the v4 region phylogeny
plate_media	one of the plate medias used to grow this organism
liquid_media	one of the liquid media used to grow this organism
notes	cultivation notes

Worksheet: full_taxonomy

This sheet contains full taxonomic information on the strains used for this manuscript. The taxonomy links it to the 'cultures' sheet. All fields are self explanatory except the following

strain_taxid	NCBI strain taxonomy ID
NCBI_Project_ID	NCBI sequencing project ID or Bioproject number
Assembly	NCBI assembly used for the construction of blast databases as described in the Supplementary Methods.

Worksheet: v4_16s

v4_16s_feature	the ID found in the 'cultures' tab
sequence	the v4 region sequence

Worksheet: full_length_16s

taxonomy	the ID found in 'full_taxonomy'
sequence	the NCBI-derived near full length 16S sequence

Worksheet: media

Medium	the name of the media used
Abbreviation	the abbreviation for the media used in the manuscript; found in Supplementary Table 5 and Supplementary Table 7
Reference	recipes or reference papers for the given media
Notes	preparation notes or modifications made to the recipes

Supplementary Table 7

Title: *in vitro* data matrices

Additional info: the worksheets in this spreadsheet are different data matrices for the *in vitro* data. These will likely be the most important data sheets for users interested in exploring the data programmatically. The rows of all the worksheets are identical – they represent a stable index that refers to the same samples across all different data types. For example, the row with index 274 is one of 3 replicates of a *Citrobacter portucalensis* (*Cp*) strain. In the 'count.raw' this row will be raw count data from *Cp* and in 'foldchange.dmrvf.fa.ps_log2' it will be log2 fold change data with media blank/delta filters, averaging of redundant molecular features, and combination of split peaks measuring the same feature. The columns of the different worksheets are similar (excluding 'aggregated_md'). They are all molecular features, or metabolites. In any worksheet with a '.fa', the columns represent metabolites; redundant molecular features have been averaged as described in the Supplementary Methods. Worksheets with a '.ps' indicate that split peaks representing the same metabolite have been summed (there are only 11 such molecular features). Finally, 'dmf' indicates a 'delta/media blank' filter, and 'dmrvf' indicates a 'delta/media blank and ratio variance' filter. Both of these are described in the Supplemental Methods.

The 'aggregated_mf' worksheet is the metadata from Supplementary Table 5 that has been aggregated so that each row refers to a physical sample measured in all three chromatographic modes. The columns descriptions are the same as those in Supplementary Table 5 and Supplementary Table 6.

The worksheets ‘regression_results’ and ‘distance_comparisons’ are the data that are used in support of Fig. 2b and Extended Data Fig. 5a. The columns are defined in detail in each worksheet.

Supplementary Table 8

Title: The *in vivo* database and *in vivo* pipeline output data matrices

Table columns:

sample_id	assigned to the physical sample
run_id	assigned to the physical sample analyzed with different methods
ms_dial_sample_name	names used in the msdial analysis output file
method columns	chromatography, ionization
sample_type	serum, urine, caecal, feces
colonization	germ-free, Bt, Cs, Bt_Ca_Er_Pd_Et, Cs_Bt_Ca_Er_Pd_Et, conventional, Cp, As
experiment	mono-colonization, community, conventional, mono-colonization_2
collection_time	denotes independent experimental repeats
mouse_id	community database-specific column, denoting individual mouse
tissue_measurement	community database-specific column, denoting three sections of the cecum
c18positive	metadata-specific column, run ids associated with c18positive method
c18negative	metadata-specific column, run ids associated with c18negative method
hilicpositive	metadata-specific column, run ids associated with hilicpositive method
metabolite columns	fc_matrix-specific: fold-change (log2-transformed) values of individual metabolite of colonized vs. germ-free controls, specific to each sample type and analytical method, one column per metabolite per analytical method, in dnames istd_corr_ion_count_matrix-specific: istd-corrected raw ion counts of individual metabolite for all sample types and colonizations states, one column per metabolite per analytical method, in dnames mode_collapsed_fc_matrix: fold-change (log2-transformed) values as described in the fc_matrix above, but with one column per metabolite. For individual metabolites detected by multiple analytical methods, their fold-change values were averaged among the preferred detection methods (See Supplementary Method for details). Column names are compound names, with co-eluting compounds separated by commas

Worksheets:

mouse sample databases: db_mono-colonization, db_community, db_conventional

in vivo sample metadata: metadata

in vivo pipeline output data matrices: raw_ion_count_matrix (columns are dnames, prior to averaging cecal section data in community experiment), istd_corr_ion_count_matrix (columns are dnames), fc_matrix (fold change, log2-transformed, in dnames), mode_collapsed_fc_matrix (fold change, log-2 transformed, columns are compound names)

Supplementary Table 9

Title: List of “phylum-associated” metabolites

Table columns:

phylum	phylum at which individual strain is associated
strain_name	bacterial strain that produced a metabolite
metabolite_produced	metabolites with ≥ 4 fold increase with corrected $P < 0.05$ when comparing their levels between bacterial supernatant vs. media controls.